

Roehl Lab Whole-Mount *in situ* Hybridization

Adapted from Thisse protocol (Nat Protoc. 2008;3(1):59-69)

Probe synthesis

Template:

set up 100ul PCR reaction with 20ng of plasmid (do not use more) and 1ul of each primer (from 25uM stock)

(Large scale:

Linearize 800ng of DNA in 20ul for 2hr.

Add 10ul each primer (100uM) to make 40ul stock

Use 1ul in 100ul PCR rxn)

96C (45sec) ---- 35X [94C (45 sec) 55C (30 sec) 72C (120sec) --- 72C (5minutes)

spin in Microcon YM-50 15min at 1000g

Add 400ul H₂O and for 45min at 1000g

Add 20ul H₂O invert and spin in new tube 1 min.

Run 1ul on a gel and quantify 1ul in spec.

Transcription mix:

Synthesis of the antisense RNA probe.

H ₂ O + 1ug DNA	21ul	
Transcription buffer (NEB)	3ul	
NTP-DIG-RNA (Boehringer)	3ul	
RNAsein (Promega)		1.5ul
T3/T7 RNA polymerase	1.5ul	

Incubate 2h at 37°C.

18ul H₂O + 2ul RNase free DNase (Roche) 30' at 37°C.

Spin sigmaspin for 15 sec at 750g

Remove base of sigmaspin and spin for 2min at 750g

Apply probe and spin for 4 min at 750g

Add 14ul RNAlater

Run 1ul on gel

Fixation and storage of embryos

(Can remove chorions by pronase treatment.)

1. Fix embryos in 4% PFA in PBS overnight to 2 days at 4°C.
2. Rinse in PBST (2 X 5') and dechorionate.

Bring through MeOH/PBST series 10' each step (30%, 60%, 3 X 100%).

4. Store in 100% MeOH at -20°C 1h (or store up to 1 year).
5. Reverse the MeOH series.
6. Rinse PBST 5 X 5'.

Digestion

1. Digest with Proteinase K (10 ug/ml) in PBST.

1 cell – 1 somite:	30 seconds
1-8 somite:	1'
9-18 somite:	3'
18 somite – 24 somite:	10'
24h (p5):	15'
30h (p13):	22'
36h/48h embryos:	30'
60h	40'
72h	50'

2. Refix in 4% PFA-PBS, 20' RT.
3. Rinse PBST 4 X 5'.

Prehyb:

1. Add 500ul Hyb+ and incubate at 70°C..
3. Incubate 2hr at 70°C.
4. Replace Hyb with 200ul of Hyb+ containing 1 to 4ul of antisense RNA probe.
5. Hybridize overnight at 70°C.

In situ Day 2:

Washes:

1. 25% 2X SSCT/75% Hyb- at 70°C 1 X 10'.
2. 50% 2X SSCT/50% Hyb- at 70°C 1 X 10'.
3. 75% 2X SSCT/25% Hyb- at 70°C 1 X 10'.
4. 2X SSCT at 70°C 1 X 10'.
5. 0.2X SSCT at 70°C 2 x 30'.
6. Allow to cool to RT 5'.
7. 25% PBST/75% 0.2X SSCT buffer at RT 10'.
8. 50% PBST/50% 0.2X SSCT buffer at RT 10'.
9. 75% PBST/25% 0.2X SSCT buffer at RT 10'.
10. Wash with PBST buffer at RT 2 X 5'.

(Optional – reflux with with 4% PFA in PBS for 20min at RT then wash well)

11. Blocking buffer 3-4h RT.
12. Incubate in antibody solution overnight with agitation at +4°C.
(anti-DIG 1:10000, 1:5000 anti-fl)

In situ Day 3:

Washes:

1. PBST at RT, brief wash.
2. PBST 6 x 15'.
3. Staining buffer 2 X 5'.

Staining:

1. Incubate embryos in staining solution at RT or 4C in the dark and monitor with a dissecting microscope.

Staining solution	1ml	3ml	5ml	15ml	25ml
NBT 50mg/ml -	4.5ul	13.5	22.5	67.5	112.5
BCIP 50 mg/ml -	3.5ul	10.5	17.5	52.5	87.5

At RT this takes 1 to 3 hrs and at 4°C leave overnight.

2. Stop the reaction by rinsing in PBST 1mM EDTA 1 X 15'.
3. Fix in 4% PFA 30'.
4. MeOH series as above and maintain in 100% MeOH 2h until blue at RT.
5. Reverse MeOH series.
6. Wash 2X in PBST 5'.
7. Bleach pigment: 30' at 37C in bleaching solution (mix every 10 min).

8. Wash 3X in PBST 5'.
9. Glycerol series (50% 10' then keep in 75-100%).
10. Store RT in the dark (up to 5 years).

Solutions:

Blocking Buffer: 50ml PBST + 1ml sheep serum + 100mg BSA.

PFA: 2gm PFA in 50ml PBS incubate at 70C 2-3 hr until dissolved (mix often) – store 1-2 days at 4C

Proteinase K (Boehringer 1000144) Make 10mg/ml suspension in PBST –20C.

20 x SSC 0.8 litre: 140.2g NaCl[3M] 65.8g NaCitrate dihydrate[0.3M] pH 7.0)

1 M Tris pH9.5 0.8litre: 95g tris base, use over 3ml HCl to pH.

4M NaCl 0.8litre: 187g NaCl

Bleaching solution: 10% H₂O₂/0.5 XSSC/5% Formamide

Heparin at 50 mg/ml (Sigma H3393) –20C.

400ml 2X Maleic buffer: 13.9g maleic acid, 4.67g NaCl, 8g NaOH pellets – check pH then adjust to then pH 7.5.

50ml staining buffer: 5ml tris pH9.5 1M, 2.5ml MgCl₂ 1M, 1.66ml NaCl 4M,, 250ul 20% Tween-20.

NBT stock: 50 mg Nitro Blue Tetrazolium in 0.7 ml of anhydrous Dimethyl-formamide + 0.3 ml H₂O 4C.

BCIP stock: 50 mg of 5-Bromo 4-Chloro3Indolyl Phosphate in 1ml anhydrous Dimethyl-formamide 4C.

Hyb+: in 50ml falcon: 25ml formamide, 25mg tRNA, 12.5ml 20X SSC, 250ul 20% Tween-20, 50ul heparin stock, 460ul 1M Citric acid and fill with H₂O.

To kill primary antibody: 3 X 5' in acidic glycine (0.1M). 50ml acidic glycine: 0.375g glycine, pH to 2.2 with concentrated HCl (about 450ul). Store at –20°C.

Fast red staining: use tris pH8.2 – one tablet in 4ml mix with p1000 in 15ml falcon, mix with 4ml of a 1:100 dilution of Naphtol stock, filter, use immediately.

Naphtol AS-MX Phosphate, Na salt: 50mg in 1ml DMSO –20C.

Holger Green: Use HRP 2° antibody and Wash in PBST. 1ml: 1.5ul fl and 2ul 3% H₂O₂ in PBST. Incubate in the dark for 30min then 4 X 10min in PBST.

Use 1% H₂O₂ (33ul/1ml) in MeOH (30min RT) to inactivate HRP antibody.

tRNA: Aliquot 800ul phenol (pH5.5) into into 18 X 2ml eppendorf tubes. Dissolve tRNA (Sigma R7876-2.5KU) in 14.5ml H₂O (RNA) and add 800ul into each of the eppendorfs. Cap and shake vigorously for 5min. Spin hard at RT for 8min. Remove supernatent almost completely to tubes containing 800ml CHCl₃. Repeat shake and spin for 4min. Remove supernatent almost completely to tubes containing 90ul 5M NH₄OAc (RNA) + 800ul isopropanol (RNA). Freeze 15 min – then spin hard at 4C for 30 min. Wash 1X with 70%EtOH (RNA) – spinning down to remove all of the wash. Dry for no more than 5min then resuspend in 2ml H₂O (RNA). Take out 1ul and mix well in 2ml H₂O and spec. Freeze in aliquots at -70C.

SIGMA R6625-25G (Ribonucleic acid from torula yeast Type VI)

or ROCHE 109 223 (Yeast RNA)

300mg in 20ml 4X SSC

10min at 68C

1X phenol/CHCl₃

Spin at 4.5K rpm in swinging bucket

1X CHCl₃

2.5ml NaOAc +25ml Isopropanol

Spin 30min

Resuspend in 4ml

Spec 1ul in 2ml (around 40ug/ml)

Aliquot 6 X 625ul (each makes 50 ml)

Black DAB - use 1% saturated $(\text{NH}_4)_2\text{Ni}_2(\text{I,II})(\text{SO}_4)_2$ in DAB mix. Store at RT.

BM purple: after washes rinse in 100mM Tris pH8.2 then BMpurple.

AB complex: Use biotinylated 2° antibody(1:200). Make AB complex 30min before use (8ul each in 1ml PBST). Incubate embryos in AB for 45min. Wash 3X in PBST.

To block endogenous AP:

- 1) after proK wash embryos briefly in H_2O , then 1hr acetic anhydride. 2.5ul acetic anhydride and 1ml 0.1M triethanolamine (pH7), mix vigorously.
- 2) use 1mM levamosol while staining

To block endogenous peroxidase:

- 1) 10min RT in 1ml 100%MeOH with 10ul 30% H_2O_2
- 2) 2% H_2O_2 in PBST 10min RT

TCA: 2% TCA is PBS 3hrs at RT, wash 3X, 4min in 0.25% Trypsin at 4C, wash 5X. 10% TCA in H_2O is stored at 4C.

ERK staining: Cells were fixed in 4% formaldehyde (fresh) for 20 min. After washes, the cells were permeabilized with 0.5% Triton X-100 in phosphate-buffered saline for 5 min, then 1hr acetic anhydride. 2.5ul acetic anhydride and 1ml 0.1M triethanolamine (pH7), mix vigorously.

1litre PBS: 1.48g NaHPO_4 , 0.43g KH_2PO_4 , 7.2 g NaCl pH to 7.2

1litre 10X PBS: 80g NaCl, 2g KCl, 14.4g $\text{Na}_2\text{HPO}_4 \cdot 7\text{H}_2\text{O}$ (9.6g $\text{Na}_2\text{HPO}_4 \cdot 2\text{H}_2\text{O}$), 2.4g KH_2PO_4 , bring up to 800ml, pH to 7.4 with HCl, then up to 1l.

400ml 2X Maleic buffer: 13.9g maleic acid, 4.67g NaCl, 8g NaOH pellets then pH 7.5. Block with 2% blocking reagent.